

# Cloning Small RNAs for Sequencing with 454 Technology

Protocol provided by Dr. Greg Hannon, Cold Spring Harbor Laboratory

## 1. RNA preparation

1. Total RNA is isolated from tissue or cells with TRIZOL® followed by two extractions with phenol chloroform, precipitation with isopropanol, and an ethanol wash.
2. RNA is then solubilized in water and the concentration is determined (should be at least 3 µg/µl)

## 2. Separation of small RNAs.

1. Prepare a 1.5 mm thick, 15% PAA urea gel with a 10-well comb. Spike radio-labelled 19mer and 24mer RNA oligos into the RNA sample, and load 50 µg of total RNA per well with 2x sample buffer in about 30-40 µl total volume. Load radiolabelled decade marker (Ambion) in flanking lanes. Run the gel until the lower dye reaches about 2/3 of gel.
2. Implant three hot gel slices and wrap the gel in Saran™ wrap. Alternatively stain the gel with SYBR Gold® (only in cases where the amount of small RNAs is sufficient to allow direct visualization). Expose the gel to a phosphoimager screen. Cut out gel slices containing the RNAs of desired size with clean scalpel. Re-expose the gel to make sure that correct sizes were excised.
3. Place the gel slices into pre-weighed siliconized Eppendorf tubes, weigh the slices and crush them with a rounded 1 ml pipette tip.
4. Add four- to six-fold volume of 0.4 M NaCl and freeze in dry ice. Thaw and shake overnight at room temperature.

## 3. Precipitation of small RNAs

1. Next day, spin the gel slices through a microcolumn (Nanosep® 100 filter) and collect the clean eluate.
2. Precipitate the small RNA by adding 1 µl of glycogen (20 µg) and 3 volumes of absolute ethanol. Put at -20°C for 4h. Spin for 30 min and remove the supernatant. Re-spin for 10 seconds and completely remove all liquid. Do not wash in 70% ethanol. Air dry the pellet for 5 min and dissolve in water (13 µl per sample).

## 4. First ligation:

1. Set up the following 20 µl ligation reaction:
  - 13 µl RNA
  - 2 µl 10X ligase buffer without ATP
  - 3 µl DMSO
  - 1 µl 3' adaptor Modban 100 µM
  - 1 µl ligase\*

\* it is best to use a mutant version of T4 RNA ligase that you prepare yourself. The reference for this protein is : Ho et al., (2004) Structure 12:327-339.

10X T4 homemade ATP-free RNA ligase buffer:

500 mM Tris-HCl pH 7.5, 100 mM MgCl<sub>2</sub>, 100 mM DTT, 600µg/ml BSA

Modban : AMP-5'p-5'p/CTGTAGGCACCATCAATdi-deoxyC- 3'  
IDT (Integrated DNA Technologies,)- Product name: miRNA  
cloning linker No. 1

2. Incubate at room temp for 1 to 2 hours.
3. Add 20 µl of 2x loading dye, heat to 65°C for a few minutes and load 40 µl of each sample per lane (10-well comb) on a 15% PAA urea gel, 1.5 mm thick. Load decade markers in flanking lanes.
4. Run the gel until the lower dye reaches about 2/3, implant hot gel slices and expose for 10 min. Excise the correct bands and elute as above, overnight.

### 5. Second ligation

1. Precipitate the ligation products as above and dissolve them each in 13 µl of water.
2. Set up the following 20 µl ligation reaction:
  - 13 µl 3' ligation product
  - 2 µl 10X T4 RNA ligase buffer (Ambion)
  - 2 µl DMSO
  - 1 µl 100 µM RNA 5' linker oligonucleotide,  
"Nelson's linker" 5'ATCGTrArGrGrCrArCrCrUrGrArArA 3'
  - 2 µl T4 RNA ligase (Ambion)

Incubate at 37°C for 1 hour.

3. Add 100 µl of 10 mM Tris (pH 8.0), 1 mM EDTA, 0.4 M NH<sub>4</sub>Acetate; Phenol extract once; Chloroform extract once.
4. Add 20 µg of glycogen; add 3 volumes (300 µl) of 100% ethanol; chill to -20°C (at least 3 hours); and centrifuge at 12,000 x g for 30 minutes at 4°C. Re-centrifuge for 10 seconds and remove all liquid. Allow the pellet to (almost) dry, and dissolve in 7 µl of water.

### 6. Reverse transcription

1. Set up the following annealing reaction:
  - 6.3 µl RNA ligation product
  - 4.2 µl 5 µM BanOne: 5'-ATTGATGGTGCCTACAG-3'
  - 
  - 10.5 µl total

Incubate at 72°C for 2 min.

Spin at 20°C for 1 min.

Cool on ice for 2 min.

2. Add the following (to make 18.9 µl total):

8.4 µl “RT Mix”	Stock of RT Mix;	30 µl 5X first strand buffer
		15 µl 20 mM DTT
		15 µl dNTPs (10 mM each)
3. Split the sample into two tubes, 9 µl in each tube:  
  
+ RT : add 1 µl (200 U) of RNase H<sup>-</sup> RT (Invitrogen Superscript™).  
- RT : add 1 µl water.
4. Incubate at 42°C for 1 hour.

### 7. PCR amplification of cDNA

1. Set up the following 100 µl PCR reaction:
  - 2 µl first-strand cDNA, or “minus RT control” (from above)
  - 74 µl water
  - 10 µl 10X PCR Buffer w/ 15 mM MgCl<sub>2</sub>
  - 10 µl dNTP Mix (2 mM each)
  - 1 µl 100 µM 5' PCR primer, BanTwo: 5'-ATCGTAGGCACCTGAAA-3'
  - 1 µl 100 µM 3' PCR primer, BanOne: 5'-ATTGATGGTGCCTACAG-3'
  - 2 µl Taq polymerase (Roche)

Mix contents by gently flicking the tube.

Centrifuge briefly to collect the contents at the bottom of the tube.

Cap the tube, and place it in a preheated (95°C) thermal cycler.

2. Run the following PCR program:
  - Step 1: 96°C 1 minute
  - Step 2: 96°C 10 sec.
  - Step 3: 50°C 60 sec
  - Step 4: 72°C 20 sec
  - Step 5: 26 cycles to Step 2
  - Step 6: 72°C 3 minutes
  - Step 7: 10°C indefinitely

### 8. Gel purification of cDNA

1. Prepare a 4% MetaPhor® Agarose gel (Cambrex) – Follow the protocol carefully for preparation of the gel.
2. Load the cDNA on the gel, alongside 4 µl of the 10 bp DNA ladder (Promega) and the rest of the “minus RT control”. Run the gel at 100 volts for 1 hour.

3. View the DNA fragments in the gel under illumination with a long wave UV source in a dark room. Mark the position of the PCR product (about 70 nt), and excise it using a clean scalpel.
4. Transfer the gel slice into a pre-weighed 1.5 ml reaction tube. If the gel slice weighs more than 250 mg, split the agarose band into two tubes. Add at least 1 volume (v/w) of 0.4 M NaCl to obtain a final volume of 500  $\mu$ l. Incubate the tube for 10 min at 70°C to melt the gel slices and add 500  $\mu$ l of 70°C pre-heated buffered water-saturated phenol (pH 7.8). Vortex the solutions vigorously and immediately separate the phases for 5 min at maximum speed in a tabletop centrifuge adjusted to room temperature. The agarose should accumulate at the interphase. Collect the aqueous upper phase and extract the aqueous phase once again with 65°C phenol.
5. For each 400  $\mu$ l of eluate, add 20  $\mu$ g of glycogen and 1 ml of ethanol. Chill at -20°C for at least one hour and recover the precipitate by centrifugation. Wash the pellet with 70% ethanol, and dissolve it in 30-50  $\mu$ l of TE buffer. *Do not allow the pellet to dry completely*, as this can cause the DNA to denature, and since it is a heterogeneous population of sequences, complementary strands may not reanneal properly. *Do not resuspend the pellet in water but in buffer for the same reason.*

Note : Alternatively the use of Qiaex II Gel Extraction Kit (Qiagen) is possible.

### 9. 454 sequencing

We start from the previous gel-purified PCR product.

1. Set up the following 100  $\mu$ l PCR reaction:
  - 1  $\mu$ l purified PCR product
  - 74  $\mu$ l water
  - 10  $\mu$ l 10X PCR Buffer w/ 15 mM MgCl<sub>2</sub>
  - 10  $\mu$ l dNTP Mix (2 mM each)
  - 1  $\mu$ l 20  $\mu$ M 5' 454 PCR primer,  
5'-GCCTCCCTCGCGCCATCAGATCGTAGGCACCTGAAA-3'
  - 1  $\mu$ l 20  $\mu$ M 3' 454 PCR primer,  
5'-GCCTTGCCAGCCCGCTCAGATTGATGGTGCCTACAG-3'
  - 2  $\mu$ l Taq polymerase (Roche)

Mix contents by gently flicking the tube.

Centrifuge briefly to collect the contents at the bottom of the tube.

2. Run the following PCR program:
  - Step 1: 96°C 1 minute
  - Step 2: 96°C 10 sec.
  - Step 3: 50°C 60 sec
  - Step 4: 72°C 20 sec
  - Step 5: 10-15 cycles to Step 2

Step 6: 72°C 3 minutes  
Step 7: 10°C indefinitely

3. Gel purify using a 2% MetaPhor® Agarose gel (or a 2% regular agarose gel) and a Qiagen gel purification kit (either Qiaex II or QIAquick)
4. Check the concentration by nanodrop A260 measurement and the quality by running 1 µl of the purified PCR product on a 2% gel. The DNA quality for 454 Sequencing has to be as good as possible.

### References:

Structure and mechanism of RNA ligase.

C. Kiong Ho, Li Kai Wang, Christopher D. Lima, and Stewart Shuman  
Structure. 2004 Feb;12(2):327-39.

Identification of novel genes coding for small expressed RNAs.

Mariana Lagos-Quintana, Reinhard Rauhut, Winfried Lendeckel, Thomas Tuschl  
Science. 2001 Oct 26;294(5543):853-8.

An Abundant Class of Tiny RNAs with Probable Regulatory Roles in *Caenorhabditis elegans*

Nelson C. Lau, Lee P. Lim, Earl G. Weinstein, David P. Bartel  
Science. 2001 Oct 26;294(5543):858-862.

An extensive class of small RNAs in *Caenorhabditis elegans*.

Rosalind C. Lee, Victor Ambros  
Science. 2001 Oct 26;294(5543):862-864.

Cloning of small RNA molecules

S. Pfeffer, M. Lagos-Quintana, T. Tuschl  
Curr. Prot. Mol. Biol, 2003, 26.4.1-26.4.18.